

# CRPS Light Therapy

Clinically tested & certified for longlasting and effective healing

## How does it work?

The regeneration of peripheral nerves, receptors and synapses is stimulated which prevents the development of neurogenic dystrophy. The volume of the focus of inflammation decreases gradually from the periphery to the centre, destroyed tissue is replaced by newly formed analogues or scar tissue. The light of CRPS Light Treatment guarantees a deep penetration in the tissues, so that the immune system gets stimulated. It reduces inflammations, stimulates blood circulation and so helps repairing tendons, ligaments and muscles.

## Why is it safe to use?

- Soft, low intensity light (2.4J)
- No known side effects
- No risk of damaging the tissues
- The light doesn't heat
- UV - free: no risk of adverse effects.

## What are the benefits?

- Reduces inflammation
- Increases micro-circulation
- Relaxes nerves and muscles
- Reduces pain transmission
- Accelerates healing regenerating cells
- Promote cell activation

## How to use it?

1. Apply the Oxygen Spray (it isn't mandatory but it cleans your skin and helps the light go deeper)
2. Place the lamp at max. 10 cm from the affected area
3. Turn the lamp on and select the suggested time (usually 10 min)
4. Repeat as many times as possible during the day (suggestion: 4 times a day minimum)